# DİREKT BİLİRUBİN

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#### REFERANSLAR

1. Rifai, N., Chiu, R. W., & Young, I., et al., (2023) Tietz Textbook of Laboratory Medicine (7th ed.), Chapter 51: Liver Disease, p.701-63.e21, Elsevier, St. Louis, Missouri 63043.
2. Pelley, J. W., PhD., (2012) Elsevier's Integrated Review Biochemistry: With Student Consult Online Access, Chapter 12: Amino Acids and Heme Metabolism, p.99-107, Elsevier Health Sciences.
3. Ferrier, D. R., (2014), Lippincott's Illustrated Reviews: Biochemistry (6th ed.), Chapter 21: Conversion of Amino Acids to Specialized Products, p.277-90, Wolters Kluwer Health.
4. Abbas, M. W., Shamshad, T., Ashraf, M. A., & Javaid, R. (2016). Jaundice: a basic review. International Journal of Research in Medical Sciences, 1313–1319. <https://dx.doi.org/10.18203/2320-6012.ijrms20161196>
5. Vendemiale G, Grattagliano I, Lupo L, Memeo V, Altomare E. Hepatic oxidative alterations in patients with extra-hepatic cholestasis. Effect of surgical drainage. J Hepatol. 2002;37(5):601-5.
6. Malhi H, Gores GJ, Malhi H, Gores GJ. Review article: the modern diagnosis and therapy of cholangiocarcinoma. Aliment Pharmacol Ther. 2006;23(9):1287-96.
7. Barkun JS, Chaudhury P, Barkun AN. Approach to the Jaundiced Patient. ACS Surgery: principles and practice. 2006.

8. Yusuf TE, Bhutani MS, Yusuf TE, Bhutani MS. Role of endoscopic ultrasonography in diseases of the extrahepatic biliary system. J Gastroenterol Hepatol. 2004;19(3):243-50.
9. Baron TH. Palliation of malignant obstructive jaundice. Gastroenterol Clin North Am. 2006;35(1):101-12.
10. Gurusamy KS, Samraj K, Gurusamy KS, Samraj K. Primary closure versus T-tube drainage after laparoscopic common bile duct stone exploration. Cochrane database of systematic reviews. 2007;(1):CD005641.
11. Gurusamy KS, Samraj K, Gurusamy KS, Samraj K. Primary closure versus T-tube drainage after open common bile duct exploration. Cochrane Database of Systematic Reviews. 2007;(1):CD005640.
12. Tai CK, Tang CN, Ha JP, Chau CH, Siu WT, Li MK. Laparoscopic exploration of common bile duct in difficult choledocholithiasis. SurgEndosc. 2004;18(6):910-4.
13. Wamsteker EJ, Wamsteker EJ. Updates in biliary endoscopy 2006. Current Opinion in Gastroenterology. 2007;23(3):324-8.
14. Strassburg, C. P. (2010). Hyperbilirubinemia syndromes (Gilbert-Meulengracht, Crigler-Najjar, Dubin-Johnson, and Rotor syndrome). Best Practice & Research in Clinical Gastroenterology, 24(5), 555–571. <https://doi.org/10.1016/j.bpg.2010.07.007>
15. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
16. VİTek, L. (2019). Bilirubin as a predictor of diseases of civilization. Is it time to establish decision limits for serum bilirubin concentrations? Archives of Biochemistry and Biophysics, 672, 108062. <https://doi.org/10.1016/j.abb.2019.108062>
17. N.N. Rehak, S.A. Cecco, G.L. Hortin, Photolysis of bilirubin in serum specimens exposed to room lighting, Clin. Chim. Acta 387 (2008) 181–183.
18. S.D. Zucker, P.S. Horn, K.E. Sherman, Serum bilirubin levels in the US population: gender effect and inverse correlation with colorectal cancer, Hepatology 40 (2004) 827–835.
19. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
20. Clinical and Laboratory Standards Institute (CLSI). Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition. CLSI Document EP34. Wayne, PA: CLSI; 2018.
21. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
22. Clinical and Laboratory Standards Institute (CLSI). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st

Edition. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.

23. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
24. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
25. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
26. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
27. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
28. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242.
29. A. Dasgupta, L.J. Langman, M. Johnson, L. Chow, Naproxen metabolites interfere with certain bilirubin assays: elimination of interference by using a Roche bilirubin assay on the Hitachi 917 analyzer, Am. J. Clin. Pathol. 133 (2010) 878–883, <https://doi.org/10.1309/AJCPN6MWATQ3SZTC>.
30. Guder WG, Narayanan S, Wisser H, et al. List of Analytes—preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag; 1996:Annex 8–9
31. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 2nd ed. Washington, DC: AACC Press; 1997:3-88.

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