# TOTAL BİLİRUBİN

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## REFERANSLAR

- Rifai, N., Chiu, R. W., & Young, I., et al., (2023) Tietz Textbook of Laboratory Medicine (7th ed.), Chapter 51: Liver Disease, p.701-763.e21, Elsevier, St. Louis, Missouri 63043.
- Pelley, J. W., PhD., (2012) Elsevier's Integrated Review Biochemistry: With Student Consult Online Access, Chapter 12: Amino Acids and Heme Metabolism, p.99-107, Elsevier Health Sciences.
- Ferrier, D. R., (2014), Lippincott's Illustrated Reviews: Biochemistry (6th ed.), Chapter 21: Conversion of

- Amino Acids to Specialized Products, p.277-90, Wolters Kluwer Health.
4. Stocker R, Yamamoto Y, McDonagh AF, Glazer AN, Ames BN. Bilirubin is an antioxidant of possible physiological importance. *Science* 1987;235:1043–1046.
  5. Mayer M. Association of serum bilirubin concentration with risk of coronary artery disease. *Clin Chem* 2000;46:1723–27.
  6. Ollinger R, Kogler P, Troppmair J, Hermann M, Wurm M, Drasche A, Königsrainer I, Amberger A, Weiss H, Ofner D, Bach FH, Margreiter R. Bilirubin inhibits tumor cell growth via activation of ERK. *Cell Cycle* 2007;6:3078–85.
  7. Schiff D, Chan G, Poznansky MJ. Bilirubin toxicity in neural cell lines N115 and NBR10A. *Pediatr Res* 1985;19:908–911.
  8. Abbas, M. W., Shamshad, T., Ashraf, M. A., & Javaid, R. (2016). Jaundice: a basic review. *International Journal of Research in Medical Sciences*, 1313–1319. <https://dx.doi.org/10.18203/2320-6012.ijrms20161196>
  9. Jacques G. Types of jaundices. *Visual Understanding Environment (VUE). Enigma*. 2009;18:55.
  10. Wickramasinghe SN, Wood WG. Advances in the understanding of the congenital dyserythropoietic anaemias. *Br J Haematol*. 2005;131:431-46.
  11. Glader B. Anemia: general consideration. In: Greer JP, eds. *Wintrobe's Clinical Hematology*. Chapter 27. Lippincott, Williams & Wilkins Co; 2004:965-75.
  12. Bektaş M, Dökmeci A, Cinar K, Halici I, Oztas E, Karayalcin S. Endoscopic management of biliary parasitic diseases. *Dig Dis Sci*. 2010;55(5):1472-8.
  13. Lidofsky SD, Kobos R. Jaundice. In: Sleisenger and Fordtran's *Gastrointestinal and Liver Disease*. 8ed. Philadelphia, Saunders Elsevier; 2006:301-16.
  14. Beckingham IJ, Ryder SD. ABC of diseases of the liver, pancreas and biliary system: investigation of liver and biliary disease. *BMJ*. 2001;322:33-6.
  15. Ryder SD, Beckingham IJ. ABC of diseases of the liver, pancreas and biliary system: other causes of parenchymal liver disease. *BMJ*. 2001;322:290-2.
  16. Roche SP, Kobos R: Jaundice in the adult patient. *Am Fam Physician*. 2004;69:299-304.
  17. Merriman RB, Peters MG. Approach to the patient with jaundice. In: Yamada T *Textbook of Gastroenterology*. 4ed. Philadelphia. Lippincott Williams and Wilkins; 2003:911-28.
  18. Mathew KG. *Medicine: Prep manual for undergraduates*. 3/e. Elsevier. India; 2008:296-7.
  19. Obeagu, E.I and Katya, M.C. (2022). A Systematic Review on Physiological Jaundice: Diagnosis and Management of the Affected Neonates. *Madonna University Journal of Medicine and Health Science*. 2 (3):25-41. <https://madonnauniversity.edu.ng/journals/index.php/medicine>
  20. Strassburg, C. P. (2010). Hyperbilirubinemia syndromes (Gilbert-Meulengracht, Crigler-Najjar, Dubin-Johnson, and Rotor syndrome). *Best Practice & Research in Clinical Gastroenterology*, 24(5), 555–571. <https://doi.org/10.1016/j.bpg.2010.07.007>
  21. Ansong-Assoku, B., D. Shah, S.D., Adnan, M., Pratibha A. Ankola, P.A.(2022). Neonatal jaundice.
  22. Vendemiale G, Grattagliano I, Lupo L, Memeo V, Altomare E. Hepatic oxidative alterations in patients with extra-hepatic cholestasis. Effect of surgical drainage. *J Hepatol*. 2002;37(5):601-5.
  23. Malhi H, Gores GJ, Malhi H, Gores GJ. Review article: the modern diagnosis and therapy of cholangiocarcinoma. *Aliment Pharmacol Ther*. 2006;23(9):1287-96.
  24. Barkun JS, Chaudhury P, Barkun AN. Approach to the Jaundiced Patient. *ACS Surgery: principles and practice*. 2006.
  25. Yusuf TE, Bhutani MS, Yusuf TE, Bhutani MS. Role of endoscopic ultrasonography in diseases of the extrahepatic biliary system. *J Gastroenterol Hepatol*. 2004;19(3):243-50.
  26. Baron TH. Palliation of malignant obstructive jaundice. *Gastroenterol Clin North Am*. 2006;35(1):101-12.
  27. Gurusamy KS, Samraj K, Gurusamy KS, Samraj K. Primary closure versus T-tube drainage after laparoscopic common bile duct stone exploration. *cochrane database of systematic reviews*. 2007;(1):CD005641.
  28. Gurusamy KS, Samraj K, Gurusamy KS, Samraj K. Primary closure versus T-tube drainage after open common bile duct exploration. *Cochrane Database of Systematic Reviews*. 2007;(1):CD005640.
  29. Tai CK, Tang CN, Ha JP, Chau CH, Siu WT, Li MK. Laparoscopic exploration of common bile duct in difficult choledocholithiasis. *SurgEndosc*. 2004;18(6):910-4.
  30. Wamsteker EJ, Wamsteker EJ. Updates in biliary endoscopy 2006. *Current Opinion in Gastroenterology*. 2007;23(3):324-8.
  31. L. Vitek, C. Bellarosa, C. Tiribelli, Induction of mild hyperbilirubinemia: hype or real therapeutic opportunity? *Clin. Pharmacol. Ther.* (2019), <https://doi.org/10.1002/cpt.1341>.
  32. F. Schiele, M. Vincent-Viry, B. Fournier, M. Starck, G. Siest, Biological effects of eleven combined oral contraceptives on serum triglycerides, gamma-glutamyltransferase, alkaline phosphatase, bilirubin and other biochemical variables, *Clin. Chem. Lab. Med.* 36 (1998) 871–878, <https://doi.org/10.1515/CCLM.1998.153>.
  33. T.W. Flaig, D.L. Gustafson, L.J. Su, J.A. Zirrolli, F. Crighton, G.S. Harrison, A.S. Pierson, R. Agarwal, L.M. Glode, A phase I and pharmacokinetic study of silybin-phytosome in prostate cancer patients, *Investig. New Drugs* 25 (2007) 139–146, <https://doi.org/10.1007/s10637-006-9019-2>.
  34. Z. Marino, G. Crespo, M. D'Amato, N. Brambilla, G. Giacobelli, L. Rovati, J. Costa, M. Navasa, X. Forns, Intravenous silybin monotherapy shows significant antiviral activity in HCV-infected patients in the peri-transplantation period, *J. Hepatol*. 58 (2013) 415–420, <https://doi.org/10.1016/j.jhep.2012.09.034>
  35. ViTek, L. (2019). Bilirubin as a predictor of diseases of civilization. Is it time to establish decision limits for

- serum bilirubin concentrations? *Archives of Biochemistry and Biophysics*, 672, 108062. <https://doi.org/10.1016/j.abb.2019.108062>
36. J. Suk, J. Jaspova, D. Biedermann, L. Petraskova, K. Valentova, V. Kren, L. Muchova, L. Vitek, Isolated Silymarin Flavonoids Increase Systemic and Hepatic Bilirubin Concentrations and Lower Lipoperoxidation in Mice, *Oxid Med Cell Long* (2019) 6026902, <https://doi.org/10.1155/2019/6026902>
  37. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1136-7
  38. Winsten S, Cehelyk B. A., rapid micro diazo technique for measuring total bilirubin. *Clin Chim Acta* 1969;25(3):441-6.
  39. Jendrassik L, Grof P. Vereinfachte photometrische Methoden zur Bestimmung des Blutbilirubins. *Biochem Z* 1938;297:81-9.
  40. Doumas BT, Poon PKC, Perry BW, et al. Candidate reference method for determination of total bilirubin in serum: development and validation. *Clin Chem* 1985;31:1779-89.
  41. Doumas BT, Perry BW, McComb RB, et al. Molar absorptivities of bilirubin (NIST SRM 916a) and its neutral and alkaline azopigments. *Clin Chem* 1990;36:1698-701.
  42. Lo SF, Jendrzyczak B, Doumas BT. Laboratory performance in neonatal bilirubin testing using commutable specimens: a progress report on a College of American Pathology study. *Arch Pathol Lab Med* 2008;132:1781-5.
  43. National Institute of Standards and Technology. Certificate of Analysis: Standard Reference Material 916a, Bilirubin. Gaithersburg, MD: National Institute of Standards and Technology; 2001.
  44. Franzini C, Cattozzo G. Low affinity complex between bilirubin and caffeine. *Clin Chem* 1987;33:597-9.
  45. Landis JB, Pardue HL. Kinetics of the reaction of unconjugated and conjugated bilirubins with p-diazobenzenesulfonic acid. *Clin Chem* 1978;24:1690-9.
  46. Lo DH, Wu TW. Assessment of the fundamental accuracy of the Jendrassik and Grof total and direct bilirubin assays. *Clin Chem* 1983;29:31-6.
  47. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
  48. N.N. Rehak, S.A. Cecco, G.L. Hortin, Photolysis of bilirubin in serum specimens exposed to room lighting. *Clin. Chim. Acta* 387 (2008) 181-183.
  49. S.D. Zucker, P.S. Horn, K.E. Sherman, Serum bilirubin levels in the US population: gender effect and inverse correlation with colorectal cancer, *Hepatology* 40 (2004) 827-835.
  50. Denney PA, Seidman DS, Stevenson DK. Drug therapy: neonatal hyperbilirubinemia. *N Engl J Med* 2001;344(8):581-90.
  51. Kaplan, M., & Hammerman, C. (2005). American Academy of Pediatrics guidelines for detecting neonatal hyperbilirubinaemia and preventing kernicterus. *Archives of Disease in Childhood-fetal and Neonatal* Rev: V1.0 Tarih: 12.2023
  - Edition, 90(6), F448-F449. <https://doi.org/10.1136/adc.2004.068726>
  52. Reddy DK, Pandey S. Kernicterus. [Updated 2023 Jun 25]. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2023 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK559120/>
  53. Bhutani, V. K., Johnson, L., & Sivieri, E. M. (1999). Predictive ability of a predischARGE hour-specific serum bilirubin for subsequent significant hyperbilirubinemia in healthy term and near-term newborns. *Pediatrics*, 103(1), 6-14. <https://doi.org/10.1542/peds.103.1.6>
  54. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
  55. Clinical and Laboratory Standards Institute (CLSI). Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition. CLSI Document EP34. Wayne, PA: CLSI; 2018.
  56. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
  57. Clinical and Laboratory Standards Institute (CLSI). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
  58. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
  59. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
  60. Bablok W et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.
  61. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
  62. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
  63. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242.
  64. A. Dasgupta, L.J. Langman, M. Johnson, L. Chow, Naproxen metabolites interfere with certain bilirubin assays: elimination of interference by using a Roche bilirubin assay on the Hitachi 917 analyzer, *Am. J. Clin. Pathol.* 133 (2010) 878-883, <https://doi.org/10.1309/AJCPN6MWWATQ3SZTC>.
  65. Quality of Diagnostic Samples, Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine, 3rd completely revised ed. 2010.

66. Data on file at Archem.

67. Jacobs DS, Oxley DK, editors. Laboratory Test Handbook, 5th ed. Hudson, OH: Lexi-Comp; 2001:117–8

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