# KOMPLEMAN C3

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## REFERANSLAR

1. Sahu, A., & Lambris, J. D. (2001). Structure and biology of complement protein C3, a connecting link between innate and acquired immunity. *Immunological Reviews*, 180, 35–48.
2. Dolmer K, Thirup S, Andersen GR, Sottrup-Jensen L, Nyborg J. Crystallization of human methylamine-treated complement C3 and C3b.
3. Einstein LP, Hansen PJ, Ballow M, et al. Biosynthesis of the third component of complement (C3) in vitro by

monocytes from both normal and homozygous C3-deficient humans. *J Clin Invest* 1977; 60: 963–9.

4. Colten HR, Ooi YM, Edelson PJ. Synthesis and secretion of complement proteins by macrophages. *Ann N Y Acad Sci* 1979; 332: 482–90.
5. Yu CL, Tsai CY, Hsieh SC, et al. Production of the third component of complement (C3) by peripheral polymorphonuclear neutrophils of the patients with rheumatoid arthritis. *Proc Natl Sci Counc Repub China B* 1995; 19: 225–32.
6. Kwan WH, van der TouwW, Heeger PS. Complement regulation of T cell immunity. *Immunol Res* 2012; 54: 247–53.
7. van Kooten C, Fiore N, Trouw LA, et al. Complement production and regulation by dendritic cells: molecular switches between tolerance and immunity. *Mol Immunol* 2008; 45: 4064–72.
8. Delanghe, J., Speeckaert, R., & Speeckaert, M. M. (2014). Complement C3 and its polymorphism: biological and clinical consequences. *Pathology*, 46(1), 1–10. <https://doi.org/10.1097/pat.0000000000000042>
9. Müller-Eberhard HJ, Schreiber RD. Molecular biology and chemistry of the alternative pathway of complement. *Adv Immunol* 1980; 29: 1–53.
10. Rodwell, V. W., Bender, D., Botham, K. M., Kennelly, P. J., & Weil, P. A. (2018b). *Harper's Illustrated Biochemistry Thirty-First Edition*. McGraw Hill Professional.
11. Rifai, N., Chiu, R. W., & Young, I., et al., (2023) *Tietz Textbook of Laboratory Medicine (7th ed.)*, Chapter 31: Amino Acids, Peptides, and Proteins, p.349-349.e42, Elsevier, St. Louis, Missouri 63043
12. Geisbrecht, B. V., Lambris, J. D., & Gros, P. (2022). Complement component C3: A structural perspective and potential therapeutic implications. *Seminars in Immunology*, 59, 101627. <https://doi.org/10.1016/j.smim.2022.101627>
13. D. Ricklin, J.D. Lambris, Complement in immune and inflammatory disorders: pathophysiological mechanisms, *J. Immunol.* 190 (2013) 3831–3838.
14. D. Ricklin, E.S. Reis, J.D. Lambris, Complement in disease: a defence system turning offensive, *Nat. Rev. Nephrol.* 2016 (12) (2016) 383–401.
15. E.S. Reis, D.C. Mastellos, D. Ricklin, A. Mantovani, J.D. Lambris, Complement in cancer: untangling an intricate relationship, *Nat. Rev. Immunol.* 18 (2018) 5–18.
16. Ram S, Lewis LA, Rice PA. Infections of people with complement deficiencies and patients who have undergone splenectomy. *Clin Microbiol Rev* 2010; 23: 740–80.
17. Tichaczek-Goska D. Deficiencies and excessive human complement system activation in disorders of multifarious etiology. *Adv Clin Exp Med* 2012; 21: 105–14.
18. Szeplaki G, Proha'szka Z, Duba J, et al. Association of high serum concentration of the third component of complement (C3) with preexisting severe coronary artery disease and new vascular events in women. *Atherosclerosis* 2004; 177: 383–9.
19. Engstrom G, Hedblad B, Eriksson KF, Janzon L, Lindgarde F. Complement C3 is a risk factor for the

- development of diabetes: a population based cohort study. *Diabetes* 2005; 54: 570–5.
20. van Oostrom AJ, Alipour A, Plokker TW, Sniderman AD, Cabezas MC. The metabolic syndrome in relation to complement component 3 and postprandial lipemia in patients from an outpatient lipid clinic and healthy volunteers. *Atherosclerosis* 2007; 190: 167–73.
  21. Wlazlo N, van Greevenbroek MM, Ferreira I, et al. Low-grade inflammation and insulin resistance independently explain substantial parts of the association between body fat and serum C3: the CODAM study. *Metabolism* 2012; 61: 1787–96.
  22. Purwar R, Langer K, Werfel T. Polymorphisms within the C3 gene are associated with specific IgE levels to common allergens and superantigens among atopic dermatitis patients. *Exp Dermatol* 2009; 18: 30–4.
  23. Tietz NW, editor. *Clinical Guide to Laboratory Tests*, 3rd ed. Philadelphia, PA: WB Saunders; 1995:164–165.
  24. Jacobs DS, DeMott WR, Grady HJ, et al. *Laboratory Test Handbook*, 4th ed. Hudson, OH: Lexi-Comp; 1996:373–374.
  25. Pickering MC, Walport MJ. Links between complement abnormalities and systemic lupus erythematosus. *Rheumatology (Oxford)* 2000; 39: 133–41.
  26. Kim, S. J., Koo, H. M., Lim, B. J., Oh, H. J., Yoo, D. E., Shin, D. H., Lee, M. J., Doh, F. M., Park, S., Yoo, T. H., Kang, S. W., Choi, K. H., Jeong, H. J., & Han, S. H. (2012). Decreased Circulating C3 Levels and Mesangial C3 Deposition Predict Renal Outcome in Patients with IgA Nephropathy. *PLOS ONE*, 7(7), e40495. <https://doi.org/10.1371/journal.pone.0040495>
  27. Augusto, J., Langs, V., Demiselle, J., Lavigne, C., Brilland, B., Duveau, A., Poli, C., Chevailler, A., Croué, A., Tollis, F., Sayegh, J., & Subra, J. (2016). Low Serum Complement C3 Levels at Diagnosis of Renal ANCA-Associated Vasculitis Is Associated with Poor Prognosis. *PLOS ONE*, 11(7), e0158871. <https://doi.org/10.1371/journal.pone.0158871>
  28. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline*. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
  29. Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
  30. Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition*. CLSI Document EP34. Wayne, PA: CLSI; 2018.
  31. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
  32. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
  33. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
  34. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
  35. Clinical and Laboratory Standards Institute (CLSI). *Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition*. CLSI Document EP37. Wayne, PA: CLSI; 2018.
  36. Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry - Third Edition*. CLSI Document EP07. Wayne, PA: CLSI; 2018.
  37. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242.
  38. Yang, P., Zhu, Z., Zang, Y., Bu, X., Xu, T., Zhong, C., Wang, A., Peng, H., Guo, D., Zheng, X., Xu, T., Chen, J., Zhang, Y., & He, J. (2021). Increased serum complement C3 levels are associated with adverse clinical outcomes after ischemic stroke. *Stroke*, 52(3), 868–877. <https://doi.org/10.1161/strokeaha.120.031715>
  39. Min, Y. G., Ju, W., Seo, J., Ha, Y., Ban, J., Kwon, Y. N., Jeong, H., Shin, J., Kim, S., Hong, Y., Kim, S. J., & Sung, J. (2023). Serum C3 complement levels predict prognosis and monitor disease activity in Guillain-Barré syndrome. *Journal of the Neurological Sciences*, 444, 120512. <https://doi.org/10.1016/j.jns.2022.120512>
  40. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002
  41. Alper CA, Rose FS. Clinical applications of complement assays. *Adv Intern Med*. 1975;20:61–88.
  42. Whicher J. Complement Component C3. In: Richie RF and Novolotskaia O, eds. *Serum Proteins in Clinical Medicine*. Scarborough: Foundation of Blood Research; 1996:1–7.

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