KOMPLEMAN C4

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REFERANSLAR

1. Tietz, N.W., Fundamentals of Clinical Chemistry, p. 940, W.B. Saunders Co., Philadelphia, 1987.
2. Tietz NW. Clinical Guide to Laboratory Test. 2nd ed. Philadelphia, PA: WB Saunders Company; 1995,52.
3. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:88-91.
4. Tietz NW, ed. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia: WB Saunders 1995:919.
5. Tietz Fundamentals of Clinical Chemistry. 5th ed. Burtis CA, Ashwood ER, eds. Philadelphia, PA: WB Saunders Company; 2001:605.
6. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
7. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
8. Clinical and Laboratory Standards Institute (CLSI). Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI Document EP17-A. Wayne, PA: CLSI; Vol. 24 No. 34.
9. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014
10. Passing-Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-79.
11. Clinical and Laboratory Standards Institute (CLSI). Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition; Approved Guideline. CLSI Document EP09-A2. Wayne, PA: CLSI; Vol. 22 No. 19.
12. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline. CLSI Document EP07. Wayne, PA: CLSI;

3rd Edition. CHERIAN G., SOLDIN ST. Clin. Chem. 27/5:748-752 (1981)

13. Young DS. Effects of Drugs on Clinical Laboratory Tests. 3rd ed. Washington: AACC Press (1990).
14. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.
15. Burtis CA, Ashwood ER, ed. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia, PA: WB Saunders Company; 1994:2183.
16. Dati F et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference range for 14 proteins in serum based on the standarization against the IFCC/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
17. Price CP, Spencer K and Whicher J. Light-scattering immunoassay of specific proteins: a review. Ann Clin Biochem 1983; 20: 1-14.
18. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.
19. Alper CA, Rose FS. Adv Intern Med. 1975;20:61-88.

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