# KREATİNİN

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### REFERANSLAR

1. Rifai, N., Chiu, R. W., & Young, I., et al., (2023) Tietz Textbook of Laboratory Medicine (7th ed.), Chapter 34: Kidney Function Tests, p.352-352.e60, Elsevier, St. Louis, Missouri 63043
2. Goldman R. Creatinine excretion in renal failure. Proc Soc Exp Biol Med 1954;85:446–8.
3. Heymsfield SB, Arteaga C, McManus C, et al. Measurement of muscle mass in humans: validity of the 24-hour urinary creatinine method. Am J Clin Nutr 1983;37:478–94.
4. Wyss M, Kaddurah-Daouk R. Creatine and creatinine metabolism. Physiol Rev 2000;80:1107–213.
5. Perrone RD, Madias NE, Levey AS. Serum creatinine as an index of renal function: new insights into old concepts. Clin Chem 1992;38:1933–53.
6. Spencer K. Analytical reviews in clinical biochemistry: the estimation of creatinine. Ann Clin Biochem 1986;23 (Pt 1):1–25.
7. Miller BF, Winkler AW. The Renal Excretion of Endogenous Creatinine in Man. Comparison with Exogenous Creatinine and Inulin. J Clin Invest 1938;17:31–40.
8. Naderer O, Nafziger AN, Bertino JS, Jr. Effects of moderatedose versus high-dose trimethoprim on serum creatinine and creatinine clearance and adverse reactions. Antimicrob Agents Chemother 1997;41:2466–70.
9. KDIGO 2012 Clinical Practice Guideline for the Evaluation and Management of Chronic Kidney Disease (Issue 1, Vol. 3). (2013). National Kidney Foundation.
10. Kidney Disease Improving Global Outcomes Acute Kidney Injury Working Group. KDIGO Clinical Practice Guideline for Acute Kidney Injury. Kidney Int 2012 suppl.; 2:1–138.

11. Jaffe M. Über den Niederschlag welchen Pikrinsäure in normalen Harn erzeugt und über eine neue Reaction des Kreatinins. Z Physiol Chem 1886;10:391-400.
12. Folin O. Beitrag zur Chemie des Kreatinins und Kreatins in Harne. Z Physiol Chem 1904;41:223-4.
13. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
14. Dilena BA. Bacterial interference with measurement of creatinine in stored plasma. Clin Chem 1988;34:1007-8.
15. Preiss DJ, Godber IM, Lamb EJ, Dalton RN, Gunn IR. The influence of a cooked-meat meal on estimated glomerular filtration rate. Ann Clin Biochem 2007;44:35-42.
16. Jacobsen FK, Christensen CK, Mogensen CE, et al. Pronounced increase in serum creatinine concentration after eating cooked meat. Br Med J 1979;1:1049–50.
17. Jacobsen FK, Christensen CK, Mogensen CE, et al. Evaluation of kidney function after meals. Lancet 1980;1:319.
18. Mayersohn M, Conrad KA, Achari R. The influence of a cooked meat meal on creatinine plasma concentration and creatinine clearance. Br J Clin Pharmacol 1983;15:227–30.
19. Pasternack A, Kuhlback B. Diurnal variations of serum and urine creatine and creatinine. Scand J Clin Lab Invest 1971;27:1–7.
20. Preiss DJ, Godber IM, Lamb EJ, et al. The influence of a cooked-meat meal on estimated glomerular filtration rate. Ann Clin Biochem 2007;44:35–42.
21. Shah KF, Stevens PE, Lamb EJ. The influence of a cooked-fish meal on estimated glomerular filtration rate. Ann Clin Biochem 2020;57:182–5.
22. Pesce, A. J., & Kaplan, L. D. (2009). Methods in Clinical Chemistry: Kaplan and Pesce's: Clinical Chemistry: Theory, Analysis, Correlation: Vol. I (5th ed.), Chapter: Creatinine, p.440-51. Elseviers.
23. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
24. Clinical and Laboratory Standards Institute (CLSI). Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition. CLSI Document EP34. Wayne, PA: CLSI; 2018.
25. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
26. Clinical and Laboratory Standards Institute (CLSI). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
27. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.

28. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
29. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
30. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
31. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
32. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242.
33. Carobene A, Ferrero C, Ceriotti F, et al. Creatinine measurement proficiency testing: assignment of matrix-adjusted ID GC-MS target values. Clin Chem 1997;43:1342-7.
34. Cobbaert CM, Baadenhuijsen H, Weykamp CW. Prime time for enzymatic creatinine methods in pediatrics. Clin Chem 2009;55:549-58.
35. Greenberg N, Roberts WL, Bachmann LM, et al. Specificity characteristics of 7 commercial creatinine measurement procedures by enzymatic and jaffe method principles. Clin Chem 2012;58:391-401
36. Weatherburn MW, Trotman RB, Jackson SH. Specific method for serum creatinine determination based on ion Exchange chromatography and an automated alkaline picrate reaction - a proposed reference method. Clin Biochem 1978;11:159-66.
37. Ali AC, Mihas CC, Campbell JA. Interferences of o-raffinose cross-linked hemoglobin in three methods for serum creatinine. Clin Chem 1997;43:1738-43.
38. Syal K, Srinivasan A, Banerjee D. Streptomycin interference in Jaffe reaction - possible false positive creatinine estimation in excessive dose exposure. Clin Biochem 2012;46:177-9.
39. Luna-Zaizar H, Virgen-Montelongo M, Cortez-Alvarez CR, et al. In vitro interference by acetaminophen, aspirin, and metamizole in serum measurements of glucose, urea, and creatinine. Clin Biochem 2015;48:538-41.
40. Swain RR, Briggs SL. Positive interference with the Jaffe reaction by cephalosporin antibiotics. Clin Chem 1977;23:1340-2.  
Kroll MH, Elin RJ. Interference with clinical laboratory analyses. Clin Chem 1994;40:1996-2005.

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