D-DİMER

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5. Adam SS, Key NS, Greenberg CS. D-dimer antigen: current concepts and future prospects. *Blood*. 2009; 113:2878–2887.
6. Linkins LA, Takach Lapner S. Review of D-dimer testing: good, bad, and ugly. *Int J Lab Hematol*. 2017;39 (Suppl 1):98–103.
7. Thachil J, Lippi G, Falavaro EJ. D-dimer testing: laboratory aspects and current issues. *Methods Mol Biol*. 2017;1646:91–104.
8. Lippi G, Cervellin G, Casagrande I, et al. D-dimer testing for suspected venous thromboembolism in the emergency department. Consensus document of AcEMC, CISMEL, SIBioC, and SIMeL. *Clin Chem Lab Med*. 2014;52:621–628.
9. Stein PD, Hull RD, Patel KC, et al. D-dimer for the exclusion of acute venous thrombosis and pulmonary embolism: a systematic review. *Ann Intern Med*. 2004;140:589–602.
10. De Monye W, Sanson BJ, Mac Gillavry MR, et al. Embolus location affects the sensitivity of a rapid quantitative D-dimer assay in the diagnosis of pulmonary embolism. *Am J Respir Crit Care Med*. 2002;165:345-348.
11. Chapman CS, Akhtar N, Campbell S, Miles K, O'Connor J, Mitchell VE. The use of D-Dimer assay by enzyme immunoassay and latex agglutination techniques in the diagnosis of deep vein thrombosis. *Clin Lab Haematol*. 1990;12:37-42.
12. Couturaud F, Kearon C, Bates SM, Ginsberg JS. Decrease in sensitivity of D-dimer for acute venous thromboembolism after starting anticoagulant therapy. *Blood Coagul Fibrinolysis*. 2002;13:241-246.
13. Reber G, De Moerloose P. Standardization of D-dimer testing. In: Kitchen S, Olson JD, Preston FE, editors. *Quality in laboratory hemostasis and thrombosis*. 2nd ed. Oxford, UK: John Wiley & Sons; 2013. p. 136–146.
14. Lip GY, Lowe GD. Fibrin D-dimer: a useful clinical marker of thrombogenesis? *Clin Sci*. 1995;89:205–214.
15. Haase C, Joergensen M, Ellervik C, et al. Age- and sex-dependent reference intervals for D-dimer: evidence for a marked increase by age. *Thromb Res*. 2013;132:676–680.
16. Nasif, W. A., Ali, A. S. E., Mukhtar, M. H., Alhuzali, A. M. H., Alnashri, Y. A., Gadah, Z. I., Edrees, E. a. A., Albarakati, H. a. M., & Aloufi, H. S. M. (2022). Elucidating the Correlation of D-Dimer Levels with COVID-19 Severity: A Scoping Review. *Anemia*, 2022, 1–9. <https://doi.org/10.1155/2022/9104209>
17. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline*. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
18. Magnette A, Chatelain M, Chatelain B, et al. Pre-analytical issues in the haemostasis laboratory: guidance for the clinical laboratories. *Thromb J*. 2016;14:49.
19. Schutgens RE, Haas FJ, Ruven HJ, et al. No influence of heparin plasma and other (pre)analytic variables on D-dimer determinations. *Clin Chem*. 2002;48: 1611–1613.
20. Toulon P, Metge S, Hangard M, et al. Impact of different storage times at room temperature of unspun citrated blood samples on routine coagulation tests results. Results of a bicenter study and review of the literature. *Int J Lab Hem*. 2017;39:458–468.
21. Bohm-Weigert M, Wissel T, Muth H, et al. Long- and short-term in vitro D-dimer stability measured with innovance D-Dimer. *Thromb Haemost*. 2010;103: 461–465.
22. Betsou F, Roussel B, Guillaume N, et al. Long-term stability of coagulation variables: protein S as a biomarker for preanalytical storage-related variations in human plasma. *Thromb Haemost*. 2009;101:1172–1175.
23. Foshat M, Bates S, Russo W, et al. Effect of freezing plasma at -20 degrees C for 2 weeks on prothrombin time, activated partial thromboplastin time, dilute Russell viper venom time, activated protein C resistance, and D-dimer levels. *Clin Appl Thromb Hemost*. 2015;21:41–47.
24. Woodhams B, Girardot O, Blanco MJ, et al. Stability of coagulation proteins in frozen plasma. *Blood Coagul Fibrinolysis*. 2001;12:229–236.
25. Wells PS, Hirsh J, Anderson DR, et al. A simple clinical model for the diagnosis of deep-vein thrombosis combined with impedance plethysmography: potential for an improvement in the diagnostic process. *J Intern Med*. 1998;243:15-23.
26. Linkins LA, Bates SM, Lang E, et al. Selective D-dimer testing for diagnosis of a first suspected episode of deep venous thrombosis: a randomized trial. *Ann Intern Med*. 2013;158:93-100.
27. Haase C, Joergensen M, Ellervik C, Joergensen MK, Bathum L. Age- and sex-dependent reference intervals for D-dimer: evidence for a marked increase by age. *Thromb Res*. 2013;132:676-680.
28. Righini M, Van Es J, Den Exter PL, et al. Age-adjusted D-dimer cutoff levels to rule out pulmonary embolism: the ADJUST-PE study. *JAMA*. 2014;311:1117-1124.
29. Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
30. Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition*. CLSI Document EP34. Wayne, PA: CLSI; 2018.
31. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
32. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
33. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
34. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative*

Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.

35. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
36. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
37. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
38. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242.
39. Mugler K, Lefkowitz JB. False-positive D-dimer result in a patient with Castleman disease. Arch Pathol Lab Med. 2004;128:328-331.
40. Huang H, Li H, Li D. Effect of serum monoclonal protein concentration on haemostasis in patients with multiple myeloma. Blood Coagul Fibrinolysis. 2015; 26:555-559.
41. Roller RE, Lahousen T, Lipp RW, et al. Elevated Ddimer results in a healthy patient. Blood Coagul Fibrinolysis. 2001;12:501-502.
42. Wu XY, Yin YF, Teng JL, et al. IgMk paraprotein from gammopathy patient can bind to cardioplipin and interfere with coagulation assay: a case report. BMC Immunol. 2017;18:32.

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