ETANOL

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REFERANSLAR

1. Tietz, N.W., Fundamentals of Clinical Chemistry, p. 940, W.B. Saunders Co., Philadelphia, 1987.
2. Tietz NW. Clinical Guide to Laboratory Test. 2nd ed. Philadelphia, PA: WB Saunders Company; 1995,52.
3. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:88-91.
4. Tietz NW, ed. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia: WB Saunders 1995:919.
5. Tietz Fundamentals of Clinical Chemistry. 5th ed. Burtis CA, Ashwood ER, eds. Philadelphia, PA: WB Saunders Company; 2001:605.
6. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
7. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
8. Clinical and Laboratory Standards Institute (CLSI). Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline.CLSI Document EP17-A. Wayne, PA: CLSI; Vol. 24 No. 34.
9. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014
10. Passing-Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-79.
11. Clinical and Laboratory Standards Institute (CLSI). Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition; Approved Guideline.CLSI Document EP09-A2. Wayne, PA: CLSI; Vol. 22 No. 19.
12. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline.CLSI Document EP07. Wayne, PA: CLSI; 3rd Edition.CHERIAN G., SOLDIN ST. Clin. Chem.

27/5:748-752 (1981)

13. Tietz Textbook of Clinical Chemistry 3.Ed. Philadelphia : W.B Saunders Company 1999 p.922-923
14. Young DS. Effects of Drugs on Clinical Laboratory Tests. 3rd ed. Washington: AACCC Press; 1990.
15. Burtis CA, Ashwood ER, eds. In: Tietz NW. Textbook of Clinical Chemistry. Toronto: WB Saunders Co; 1994:1170-1175.
16. Thomas L. Clinical Laboratory Diagnostics. 1. Ed. Frankfurt: TH-Books Verlagsgesellschaft ; 1998 p. 1168-1170
17. Lundquist F. The Determination of Ethyl Alcohol in Blood and Tissues. In: Glick D, ed. Methods Biochem Anal. Vol III. New York, NY: Interscience publishers; 1957:217- 251.
18. William H. Porter Ph.D.Clinical Toxicology In: Burtis CA, Ashwood ER, editorts.
19. Bonnicksen RK, Theorell H. An Enzymatic Method for the Micro Determination of Ethanol. Scand J Clin Lab Invest. 1951;3:58.
20. Clinical and Laboratory Standards Institute (formerly NCCLS). How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory; Approved Guideline. Wayne, PA: NCCLS; 1995. Clinical and Laboratory Standards Institute Document C29-A.
21. Clinical and Laboratory Standards Institute (formerly NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. NCCLS Document EP05-A2.

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