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4. Chang, J.; Bird, R.; Clague, A.; Carter, A. Clinical utility of serum soluble transferrin receptor levels and comparison with bone marrow iron stores as an index for iron deficient erythropoiesis in a heterogeneous group of patients. *Pathology* 2007, 39, 349–353.
5. Wish, J.B. Assessing iron status: Beyond serum ferritin and transferrin saturation. *Clin. J. Am. Soc. Nephrol.* 2006, 1, S4–S8.
6. Aisen P, Leibman A, Reich HA. Studies on the binding of iron to transferrin and conalbumin. *J Biol Chem* 1966;241:1666–71.
7. Gambino, R.; Desvarieux, E.; Orth, M.; Matan, H.; Ackattupathil, T.; Lijoi, E.; Wimmer, C.; Bower, J.; Gunter, E. The relation between chemically measured total iron-binding capacity concentrations and immunologically measured transferrin concentrations in human serum. *Clin. Chem.* 1997, 43, 2408–2412.
8. Guo, R., Gao, J., Hui, L., et al., (2022). An Improved Method for Quick Quantification of Unsaturated Transferrin. *Biosensors*, 12(9), 708. <https://doi.org/10.3390/bios12090708>
9. Singh, M., Patra, S., & Singh, R. (2021). Common techniques and methods for screening of natural products for developing of anticancer drugs. In Elsevier eBooks (pp. 323–353). <https://doi.org/10.1016/b978-0-12-821710-8.00015-1>
10. Persijn, J.P.; van der Slik, W.; Riethorst, A. Determination of serum iron and latent iron-binding capacity (LIBC). *Clin. Chim. Acta* 1971, 35, 91–98.
11. Faruqi A, Mukkamalla SKR. Iron Binding Capacity. 2023 Jan 2. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2023 Jan–. PMID: 32644545.
12. Pfeiffer CM, Looker AC. Laboratory methodologies for indicators of iron status: strengths, limitations, and analytical challenges. *Am J Clin Nutr* 2017;106:1606S–14S.
13. Ramsay WN. The determination of the total iron-binding capacity of serum. *Clin Chim Acta* 1957;2:221–6.
14. The measurement of total and unsaturated iron-binding capacity in serum. *Br J Haematol* 1978;38:281–7.
15. Revised recommendations for the measurements of the serum iron in human blood. Iron Panel of the International Committee for Standardization in Haematology. *Br J Haematol* 1990;75:615–6.
16. Yamanishi H, Kimura S, Iyama S, Yamaguchi Y, Yanagihara T. Fully automated measurement of total iron-binding capacity in serum. *Clin Chem* 1997;43:2413–7.
17. Siek G, Lawlor J, Pelczar D, Sane M, Musto J. Direct serum total iron-binding capacity assay suitable for automated analyzers. *Clin Chem* 2002;48:161–6.
18. Gambino R. Serum transferrin (total iron binding capacity) in evaluation of iron status. *Clin Chem* 1996;42:2053.
19. Pesce, A. J., & Kaplan, L. D. (2009). *Methods in Clinical Chemistry: Kaplan and Pesce's: Clinical Chemistry: Theory, Analysis, Correlation: Vol. I* (5th ed.), Chapter: Iron and Iron Binding Capacity, p.781-84. Elseviers.
20. Åsberg, A., Thorstensen, K., Mikkelsen, G., & Åsberg, A. (2013). The diagnostic accuracy of unbound iron binding capacity (UIBC) as a test for empty iron stores. *Scandinavian Journal of Clinical & Laboratory Investigation*, 73(3), 208–213. <https://doi.org/10.3109/00365513.2013.765029>
21. Koseoglu, M.; Hur, A.; Atay, A.; Cuhadar, S. Effects of hemolysis interferences on routine biochemistry parameters. *Biochem. Med.* 2011, 21, 79–85
22. Shin, D.H.; Kim, J.; Uh, Y.; Lee, S.I.; Seo, D.M.; Kim, K.S.; Jang, J.Y.; Lee, M.H.; Yoon, K.R.; Yoon, K.J. Development of an integrated reporting system for verifying hemolysis, icterus, and lipemia in clinical chemistry results. *Ann. Lab. Med.* 2014, 34, 307–312.
23. Tavill, A. S. (2001). Diagnosis and management of hemochromatosis. *Hepatology*, 33(5), 1321–1328. <https://doi.org/10.1053/jhep.2001.24783>
24. Rayman, M. P., Barlis, J., Evans, R., Redman, C. W., & King, L. J. (2002). Abnormal iron parameters in the pregnancy syndrome preeclampsia. *American Journal of Obstetrics and Gynecology*, 187(2), 412–418. <https://doi.org/10.1067/mob.2002.123895>
25. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
26. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
27. Clinical and Laboratory Standards Institute (CLSI). Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition. CLSI Document EP34. Wayne, PA: CLSI; 2018
28. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
29. Clinical and Laboratory Standards Institute (CLSI). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
30. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
31. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
32. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in

33. Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
34. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
35. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242).
36. Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders 1987:789-824.

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