DEMİR

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REFERANSLAR

1. Rifai, N., Chiu, R. W., & Young, I., et al., (2023) Tietz Textbook of Laboratory Medicine (7th ed.), Chapter 40: Iron Metabolism, p.418-e40, Elsevier, St. Louis, Missouri 63043
2. Pantopoulos K, Porwal SK, Tartakoff A, Devireddy L. Mechanisms of mammalian iron homeostasis. *Biochemistry* 2012;51:5705–24.
3. Koppenol WH. The centennial of the Fenton reaction. *Free Radic Biol Med* 1993;15:645–51.
4. Torti FM, Torti SV. Regulation of ferritin genes and protein. *Blood* 2002;99:3505–16.
5. Drakesmith H, Prentice AM. Hcpidin and the iron-infection axis. *Science* 2012;338:768–72.
6. Ganz T. Hcpidin and iron regulation, 10 years later. *Blood* 2011;117:4425–33.
7. Nemeth E, Tuttle MS, Powelson J, Vaughn MB, Donovan A, Ward DM, Ganz T, et al. Hcpidin regulates cellular iron efflux by binding to ferroportin and inducing its internalization. *Science* 2004;306:2090–3.
8. Ganz T. Systemic iron homeostasis. *Physiol Rev* 2013;93:1721–41.
9. Kuhn LC. Iron regulatory proteins and their role in controlling iron metabolism. *Metallomics* 2015;7:232–43.
10. Finch C. Regulators of iron balance in humans. *Blood* 1994;84:1697–702.
11. Cook JD, Skikne BS, Lynch SR, Reusser ME. Estimates of iron sufficiency in the US population. *Blood* 1986;68:726–31.
12. Arosio P, Ingrassia R, Cavadini P. Ferritins: a family of molecules for iron storage, antioxidation and more. *Biochim Biophys Acta* 2009;1790:589–99.
13. Kassebaum NJ, Jasrasaria R, Naghavi M, Wulf SK, Johns N, Lozano R, Regan M, et al. A systematic analysis of global anemia burden from 1990 to 2010. *Blood* 2014;123:615–24.
14. Stevens GA, Finucane MM, De-Regil LM, Paciorek CJ, Flaxman SR, Branca F, Pena-Rosas JP, et al. Global, regional, and national trends in haemoglobin concentration and prevalence of total and severe anaemia in children and pregnant and non-pregnant women for 1995–2011: a systematic analysis of populationrepresentative data. *Lancet Glob Health* 2013;1:e16–25.
15. WHO, UNICEF, UNU. Iron deficiency anaemia: assessment, prevention, and control. A guide for programme managers. Geneva, World Health Organisation 2001;WHO/NHD/01.3.

- Available at http://www.who.int/nutrition/publications/micronutrients/anaemia_iron_deficiency/WHO_NHD_01.3/en/.
16. Camaschella C. Iron-deficiency anemia. *N Engl J Med* 2015;372:1832–43.
 17. +Pasricha SR, Drakesmith H, Black J, Hipgrave D, Biggs BA. Control of iron deficiency anemia in low- and middle-income countries. *Blood* 2013;121:2607–17.
 18. Bothwell TH. Iron requirements in pregnancy and strategies to meet them. *Am J Clin Nutr* 2000;72:257S–64S.
 19. Bruner AB, Joffe A, Duggan AK, Casella JF, Brandt J. Randomised study of cognitive effects of iron supplementation in non-anaemic iron-deficient adolescent girls. *Lancet* 1996;348:992–6.
 20. Beard JL, Hendricks MK, Perez EM, Murray-Kolb LE, Berg A, Vernon-Feagans L, Irlam J, et al. Maternal iron deficiency anemia affects postpartum emotions and cognition. *J Nutr* 2005;135:267–72.
 21. Beard JL. Iron biology in immune function, muscle metabolism and neuronal functioning. *J Nutr* 2001;131:568S–79S; discussion 580S.
 22. Ahluwalia N, Sun J, Krause D, Mastro A, Handte G. Immune function is impaired in iron-deficient, homebound, older women. *Am J Clin Nutr* 2004;79:516–21.
 23. DeLoughery TG. Microcytic anemia. *N Engl J Med* 2014;371:2537.
 24. Pratt JJ, Khan KS. Non-anaemic iron deficiency - a disease looking for recognition of diagnosis: a systematic review. *Eur J Haematol* 2015.
 25. Pasricha SR, Low M, Thompson J, Farrell A, De-Regil LM. Iron supplementation benefits physical performance in women of reproductive age: a systematic review and meta-analysis. *J Nutr* 2014;144:906–14.
 26. Avni T, Reich S, Lev N, Gafter-Gvili A. Iron supplementation for restless legs syndrome - A systematic review and metaanalysis. *Eur J Intern Med* 2019;63:34–41.
 27. Pasricha SR, Tye-Din J, Muckenthaler MU, Swinkels DW. Iron deficiency. *Lancet* 2021;16:233–48.
 28. Low M, Farrell A, Biggs BA, Pasricha SR. Effects of daily iron supplementation in primary-school-aged children: systematic review and meta-analysis of randomized controlled trials. *CMAJ* 2013;185:E791–802.
 29. Klip IT, Comin-Colet J, Voors AA, Ponikowski P, Enjuanes C, Banasiak W, Lok DJ, et al. Iron deficiency in chronic heart failure: an international pooled analysis. *Am Heart J* 2013;165:575–82 e573.
 30. Thomas DW, Hinchliffe RF, Briggs C, Macdougall IC, Littlewood T, Cavill I, British Committee for Standards in H. Guideline for the laboratory diagnosis of functional iron deficiency. *Br J Haematol* 2013;161:639–48.
 31. Weiss G, Goodnough LT. Anemia of chronic disease. *N Engl J Med* 2005;352:1011–23.
 32. Gasche C, Berstad A, Befrits R, Beglinger C, Dignass A, Erichsen K, Gomollon F, et al. Guidelines on the diagnosis and management of iron deficiency and anemia in inflammatory bowel diseases. *Inflamm Bowel Dis* 2007;13:1545–53.
 33. Weiss G, Ganz T, Goodnough LT. Anemia of inflammation. *Blood* 2019;133:40–50.
 34. Fleming RE, Ponka P. Iron overload in human disease. *N Engl J Med* 2012;366:348–59
 35. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
 36. Pilon VA, Howanitz PJ, Howanitz JH, Domres N. Day-to-day variation in serum ferritin concentration in healthy subjects. *Clin Chem* 1981;27:78–82.
 37. Borel MJ, Smith SM, Derr J, Beard JL. Day-to-day variation in iron-status indices in healthy men and women. *Am J Clin Nutr* 1991;54:729–35.
 38. Maes M, Bosmans E, Scharpe S, Hendriks D, Cooremans W, Neels H, De Meyer F, et al. Components of biological variation in serum soluble transferrin receptor: relationships to serum iron, transferrin and ferritin concentrations, and immune and haematological variables. *Scand J Clin Lab Invest* 1997;57:31–41.
 39. Statland BE, Winkel P. Relationship of day-to-day variation of serum iron concentrations to iron-binding capacity in healthy young women. *Am J Clin Pathol* 1977;67:84–90.
 40. Dallman PR. Diagnosis of anemia and iron deficiency: analytic and biological variations of laboratory tests. *Am J Clin Nutr* 1984;39:937–41.
 41. Hershko C, Bar-Or D, Gaziel Y, Naparstek E, Konijn AM, Grossowicz N, Kaufman N, et al. Diagnosis of iron deficiency anemia in a rural population of children. Relative usefulness of serum ferritin, red cell protoporphyrin, red cell indices, and transferrin saturation determinations. *Am J Clin Nutr* 1981;34:1600–10.
 42. Ridefelt P, Larsson A, Rehman JU, Axelsson J. Influences of sleep and the circadian rhythm on iron-status indices. *Clin Biochem* 2010;43:1323–8.
 43. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
 44. Clinical and Laboratory Standards Institute (CLSI). Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking – 1st Edition. CLSI Document EP34. Wayne, PA: CLSI; 2018.
 45. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
 46. Clinical and Laboratory Standards Institute (CLSI). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
 47. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
 48. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.

49. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
50. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
51. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
52. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131:41194-242.
53. Wu AHB, editor. Tietz Clinical Guide to LaboratoryWuy Tests. 4th ed. St. Louis, MO: Elsevier Saunders; 2006
54. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
55. Seligman PA, Schleicher RB. Comparison of methods used to measure serum iron in the presence of iron gluconate or iron dextran. Clin Chem 1999;45:898–901.
56. Huisman W. Interference of imferon in colorimetric assays for iron. Clin Chem 1980;26:635–7.
57. Ikuta K, Ito S, Tanaka H, Sasaki K, Torimoto Y, Fujiya M, Kohgo Y. Interference of deferasirox with assays for serum iron and serum unsaturated iron binding capacity during iron chelating therapy. Clin Chim Acta 2011;412:2261–6.
58. Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.

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SYMBOLS

IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
R1	Reagent 1
R2	Reagent 2
GTIN	Global Trade Item Number
REF	Reference Number
GLP	Good Laboratory Practices
FOR USE WITH	Identifies Products to Be Used Together
PRODUCT OF TURKEY	Product of Turkey
	Manufacturer
	Expiration Date
	Temperature Limits
	Consult Instructions for Use
	Caution
	Number of Tests